

Biofiller Plasma Gel

Model # :AGM-10

- Add volume and smooth out the deep creases and wrinkles in the face
- Plumping thin lips, reduce fine lines around the neck
- Augment cheeks to enhance youthfulness, reduce or remove eye-bags
- Improve the appearance of recessed scars, enhance shallow contours
- Turn Platelet-Poor Plasma (PPP) into Plasma Bio-Filler Gel
- Injected on the face, neck, breasts and butt to correct lack of volume, wrinkles, deep folds, unaesthetic lines, and unusual thinness.



Technical Specification

Model		AGM-10
Temperature range		25-90°C
Max Volume		10 Syringes
Timer		0-99min
Record function		Yes
Voltage		220V 50Hz / 110V 60Hz
Dimension		290×280×140mm
Weight		6.5kg
Cat No.	110V 60HZ	70991100
	220V 60HZ	70992100